

ENDOWRIST® 8 MM LONG BIPOLAR GRASPER FOR THE DA VINCI XI® SURGICAL SYSTEM



Optimized for **ERBE**
AutoStop feature

Blunt tip for dissection

Improved tissue
grasping*

WHAT SURGEONS ARE SAYING

“You got it totally right with this one. It picks up tissue really well, doesn't tear, there's minimal thermal spread, and no slipping.†”

Dr. Michael Stifelman | URO

“This instrument is an all-rounder. It coagulates, blunt dissects, and has a wide jaw open angle. I have a new bipolar now-the Long Bipolar Grasper!†”

Dr. Kathy Huang | GYN

ORDERING INFORMATION







PN 470400

10 uses

*As compared to *da Vinci Si*® PK Dissecting Forceps

†Comments reflect an individual surgeon experience that may or may not be reproducible and is not generalizable.

NEW *ENDOWRIST*® 8 MM LONG BIPOLAR GRASPER

	<i>Xi</i>	<i>Si</i>
	LONG BIPOLAR GRASPER	MARYLAND DISSECTING FORCEPS
JAW OPEN ANGLE	 70°	 40°
BLUNT DISSECTION	 1.7 CM	 1.5 CM
		PK DISSECTING FORCEPS
		 70°
		 1.6 CM

INTUITIVE
SURGICAL®

IMPORTANT SAFETY INFORMATION

Serious complications may occur in any surgery, including *da Vinci*® Surgery, up to and including death. Examples of serious or life-threatening complications, which may require prolonged and/or unexpected hospitalization and/or reoperation, include but are not limited to one or more of the following: injury to tissues/organs, bleeding, infection and internal scarring that can cause long-lasting dysfunction/pain. Individual surgical results may vary.

Risks specific to minimally invasive surgery, including *da Vinci*® Surgery, include but are not limited to, one or more of the following: temporary pain/nerve injury associated with positioning; a longer operative time, the need to convert to an open approach, or the need for additional or larger incision sites. Converting the procedure could result in a longer operative time, a longer time under anesthesia, and could lead to increased complications. Contraindications applicable to the use of conventional endoscopic instruments also apply to the use of all *da Vinci* instruments. You should discuss your surgical experience and review these and all risks with your patients, including the potential for human error and equipment failure. Physicians should review all available information. Clinical studies are available through the National Library of Medicine at www.ncbi.nlm.nih.gov/pubmed.

Be sure to read and understand all information in the applicable user manuals, including full cautions and warnings, before using *da Vinci* products. Failure to properly follow all instructions may lead to injury and result in improper functioning of the device. Training provided by Intuitive Surgical is limited to the use of its products and does not replace the necessary medical training and experience required to perform surgery. Procedure descriptions are developed with, reviewed and approved by independent surgeons. Other surgical techniques may be documented in publications available at the National Library of Medicine. For Important Safety Information, indications for use, risks, full cautions and warnings, please also refer to www.davincisurgery.com/safety and www.intuitivesurgical.com/safety. Unless otherwise noted, products featured are available for commercial distribution in the U.S. For availability outside the U.S., please check with your local representative or distributor.

© 2016 Intuitive Surgical, Inc. All rights reserved. Product names are trademarks or registered trademarks of their respective holders.
PN 1030430-US Rev A 11/16